

SUMMER OFFERS

01.05-31.08 2024

Don't miss out on these limited-time deals!



Upgrade your lab today.

Do you need to upgrade your lab setup? You can now find an exclusive selection of basic instruments at unbeatable prices this spring and summer from our Dominique Dutscher portfolio!

490005 Plexiglas DNA and PCR workstation

Ideal for safe sample preparation. UV radiation disinfects the cabinet work surface, inactivating any DNA/RNA within 15-30 min / Radiation automatically switched off when the door is opened. Bactericidal cabinet decontamination system, by air cleaning-recirculation and ozone-free radiation.

List price: ~~24.161 DKK~~

Campaign price: 17.187 DKK



486074 PCR-3A PCR cabinet

Optimal safety: Complies with norm ISO 14644-1 class 4. Vertical laminar flow of clean air to guarantee sample protection. Pre-filtered beforehand to optimise the life span of the filter. High power UV tube for destruction of DNA or RNA fragments.

List price: ~~51.298 DKK~~

Campaign price: 30.218 DKK



Airstream E series Microbiological Safety Station

Workstation with 2 electric plugs, HEPA filter and a UV lamp and antimicrobial ISOCIDE coating. CE/EN compliant Certified for AC2-3E8, AC2-4E8, AC2-6E8.

- Comes in 90 cm, 120 cm, and 180 cm width
- Stand sold separately



Campaign price:

486343: 90 cm	486344: 120 cm	486346: 180 cm
71.891 DKK	73.946 DKK	101.532 DKK
47.921 DKK	47.921 DKK	66.152 DKK

When ordering from this campaign, please use campaign code: CM00431

082137B Arpege 140 storage dewar

Medium capacity container for the cryopreservation of small biological products. Useful capacity (l) : 144, Capacity (with racks) 4.800 2 ml tubes.

List price: 52.369 DKK

Campaign price: 39.037 DKK



82733 Arpege liquid air 140 cryopreserver - optimal pack 3

Medium capacity container for the cryopreservation of small biological products with CryoMemo, automatic filling and degassing. Intuitive menu interface, allowing access to all data - level and temperature display, alarm settings, etc. via Internet server for remote monitoring.

List price: 129.544 DKK

Campaign price: 100.052 DKK



333352 Ohaus Pioneer analytical balance PX224

Balance offering accuracy and repeatability in essential weighing applications in laboratory. RS232 port for data transfer. USB port for data extraction and GLP printing option (date, time, identifiers).

List price: 14.687 DKK

Campaign price: 12.470 DKK



092689 Ohaus Explorer semi-micro balance EX125

Balance with large VGA color touch-screen. Base and display are separated for flexibility and display can be mounted on a tower or wall. GLP-compliant printouts (date, hour, identification) Precision 0.01 mg.

List price: 56.700 DKK

Campaign price: 48.141 DKK



144783 OLS 26 - shaking water bath

Combined orbital/linear shaking water bath with precision digital temperature control. Includes universal tray and polycarbonate lid as well as drain tap.

List price: 56.968 DKK

Campaign price: 35.118 DKK



158310 Integra Vacusafe Comfort

For completely safe suction without contamination of supernatant liquids in a 4 litre polypropylene flask with level detector and -300 to -600 mbar vacuum.

All accessories (handle, connectors, tube, bottle) are autoclavable.

Includes all elements of Vacusafe Comfort ECO plus lid with rapid connection and level detector.

List price: 17.577 DKK

Campaign price: 14.140 DKK



255201 TurboCycler lite Gradient thermocycler

96 well gradient block compatible with PCR 0.2 ml standard tubes and «low profile» tubes and strips, skirted, semi-skirted and non-skirted plates. Monochrome screen and intuitive interface via touch button control.

List price: 33.681 DKK

Campaign price: 25.062 DKK



255200 TurboCycler 2 Gradient thermocycler

96 well gradient block compatible with standard and «low profile» tubes and strips, and skirted, semi-skirted and non-skirted plates. Intuitive control interface with icons and graphics via 7" color touch screen that can be manipulated with laboratory gloves.

List price: 53.857 DKK

Campaign price: 34.031 DKK



TCTT-4810 TurboCycler 3 Thermal Cycler

Run 3 independent experiments (3x48 samples) simultaneously with unmatched flexibility. Fast gradient optimization and advanced ramp rate control tailor PCR conditions for maximized efficiency.

List price: 75.000 DKK

Campaign price: 56.250 DKK



531623 Magnetapure mini nucleic acid extraction automat

Open system for 16 parallel magnetic bead purifications DNA, RNA, proteins and cells from a wide variety of samples such as blood, bacteria, tissue, body fluids and plants. Easy programming to adapt your current magnetic beads or contact us to learn already validated protocols.

List price: 58.984 DKK

Campaign price: 41.869 DKK



255744 Magnetapure 32 Plus nucleic acid extraction automat

Open system for 32 parallel magnetic bead purifications DNA, RNA, proteins and cells from a wide variety of samples such as blood, bacteria, tissue, body fluids and plants. Easy programming to adapt your current magnetic beads or contact us to learn already validated protocols.

List price: 141.750 DKK

Campaign price: 85.260 DKK



255745 Magnetapure 96 Plus nucleic acid extraction automat

Open system for 96 parallel magnetic bead purifications DNA, RNA, proteins and cells from a wide variety of samples such as blood, bacteria, tissue, body fluids and plants. Easy programming to adapt your current magnetic beads or contact us to learn already validated protocols.

List price: 270.000 DKK


Campaign price: 175.691 DKK




Looking for something else?


Our selection of products span the whole range of needs for your lab - from essential consumables to large research- and diagnostics instruments - and everything in between.


Browse our web shop or even easier; contact your local Customer Care team and they can help you.



SCAN or CLICK

 +45 8745 9010

 order@ahdiagnostics.dk

 www.ahdiagnostics.dk

When ordering from this campaign, please use campaign code: CM00431

