

END-OF-YEAR OFFERS

FROM EPPENDORF

Save 15 % on essential lab equipment from Eppendorf

ThermoMixer® C Exceeds Your Expectations

Mixing, heating, and cooling in one instrument for high flexibility. Accommodates all common vessels and plate formats from 5 µL to 50 mL using Eppendorf Smartblocks™. For a condensation-free sample, use the Eppendorf ThermoTop®.

Used in various applications such as:

- Plasmid/RNA/DNA purification
- cDNA synthesis
- Extraction of DNA fragments from agarose gels
- Enzymatic reactions
- Bacterial growth in microreaction vessels
- Labeling of DNA/RNA and proteins
- Lysis reaction at 100C



CellXpert® C170i ensures safe growth conditions - even for your sensitive cells

- Optimize your growth conditions with a CO₂-incubator that have a fast temperature and CO₂ recovery in less than 5 min, a fanless design and good temperature homogeneity
- Easy cleaning and reliable contamination prevention by 180 °C sterilization in a seamless chamber with the option for 4- or 8-segmented inner doors or an antimicrobial copper interior
- Precise documentation and track of the CO₂, temperature, door opening times ect. for easy troubleshooting in the case of an alarm notification.



Offers are valid from 01.09.2024 until 31.10.2024
Please use campaign code when ordering: CM00443

Eppendorf Research® Plus Pipette - an Essential Tool in Your Lab

The ultra-light **Eppendorf Research Plus** manual pipette not only meets the highest needs in precision and accuracy, but it also provides you with ultimate ergonomics and increased flexibility.

- 3123 000.900 3-pack, 0.5-10 µL/10-100 µL/100-1000 µL, incl. epT.I.P.S.® boxes
3123 000.918 3-pack, 2-20 µL/20-200µL/100-1000 µL, incl. epT.I.P.S.® boxes
3123 000.926 3-pack, 100-1000 µL/0,5-5 mL/1-10 mL, incl. epT.I.P.S.® boxes



Research® plus Move It® Easy Adjustment of your pipetting

- Multi-channel pipette with adjustable tip spacing and a 360° rotating pipette head
- Quick and easy changeover of different plate formats thanks to predefined parameters
- Suitable for use with plates (384, 96, 48, 24 and 12 wells), tubes (1.5 and 2 ml) and agarose gels
- Optimised ergonomics and balance to reduce muscle fatigue
- Easily autoclavable




Offers are valid from 01.10.2024 until 15.11.2024
Please use campaign code when ordering: CM00444

Looking for something else?

Our selection of products span the whole range of needs for your lab - from essential consumables to large research- and diagnostics instruments - and everything in between.

Browse our web shop or even easier; contact your local Customer Care team and they can help you.



+45 8745 9010



order@ahdiagnostics.dk



www.ahdiagnostics.dk

